

## *General Protocol & Storage*

### Product Description

A set of Stellaris RNA FISH Probes is comprised of up to 48 singly labeled oligonucleotides designed to selectively bind to targeted transcripts. Stellaris RNA FISH Probes bound to target RNA produce fluorescent signals that permit detection of single RNA molecules as diffraction-limited spots by conventional fluorescence microscopy.

### Storage Guidelines

#### *Stellaris RNA FISH Probes*

Stellaris RNA FISH Probes are shipped dry and can be stored at +2 to +8 °C in this state. Dissolved probe mix should be subjected to a minimum number of freeze-thaw cycles. For daily and short-term use of dissolved probe mix, storage at +2 to +8 °C in the dark for up to a month is recommended. For storage lasting longer than a month, we recommend aliquoting and freezing probes in the dark at -15 to -30 °C.

#### *Stellaris RNA FISH Hybridization Buffer*

Stellaris RNA FISH Hybridization Buffer should be stored at +2 to +8 °C for short-term and long-term use.

#### *Stellaris RNA FISH Wash Buffer A and Wash Buffer B*

Stellaris RNA FISH Wash Buffers A and B should be stored at room temperature for short-term and long-term use.

### Reagents and Equipment

- a) TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0)
- b) 37% Formaldehyde Solution
- c) 10X Phosphate Buffered Saline (PBS), RNase-free
- d) Nuclease-free water
- e) Deionized Formamide
- f) Ethanol for molecular biology
- g) Stellaris RNA FISH Hybridization Buffer (Biosearch Technologies Cat# SMF-HB1-10)
- h) Stellaris RNA FISH Wash Buffer A (Biosearch Technologies Cat# SMF-WA1-60)
- i) Stellaris RNA FISH Wash Buffer B (Biosearch Technologies Cat# SMF-WB1-20)
- j) 4',6-diamidino-2-phenylindole (DAPI)
- k) Vectashield<sup>®</sup> Mounting Medium (Vector Laboratories Cat #H-1000)
- l) 2-well chambered coverglass
- m) 18 x 18 mm square #1 coverglass
- n) RNase free consumables such as pipette tips
- o) Humidified chamber (or equivalent): 150 mm tissue culture plate; place a single water-saturated paper towel alongside the inner chamber edge
- p) 37 °C laboratory oven

### Microscope:

- Wide-field fluorescence microscope (e.g., Nikon Eclipse Ti or equivalent). We provide limited support for confocal applications.
- A high numerical aperture (>1.3) and 60-100x oil-immersion objective.
- Strong light source, such as a mercury or metal-halide lamp (newer LED-based light sources may also be sufficient).
- Filter sets appropriate for the fluorophores.
- Standard cooled CCD camera, ideally optimized for low-light level imaging rather than speed (13 µm pixel size or less is ideal).

### Preparation of Reagents

**NOTE:** When performing Stellaris RNA FISH, it is imperative to limit RNA degradation. Please ensure that all consumables and reagents are RNase-free. Recipes below are for set volumes. Please adjust accordingly.

#### Reconstituting the dried probe stock:

**ShipReady Probe Set (1 nmol):** A ShipReady probe set can provide up to 80 hybridizations. Re-dissolve the dried oligonucleotide probe blend in 80 µL of TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0) to create a probe stock of 12.5 µM. *Mix well by pipetting up and down*, and then vortex and centrifuge briefly.

**DesignReady or Custom Probe Set (5 nmol):** A DesignReady or custom probe set can provide up to 400 hybridizations. Re-dissolve the dried oligonucleotide probe blend in 400 µL of TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0) to create a probe stock of 12.5 µM. *Mix well by pipetting up and down*, and then vortex and centrifuge briefly.

#### Fixation Buffer:

*Final composition is 3.7% (vol./vol.) formaldehyde in 1X PBS*

**For a final volume of 10 mL, mix:**

- 1 mL 37% Formaldehyde solution
- 1 mL 10X Phosphate Buffered Saline (PBS), RNase-free
- 8 mL Nuclease-free water

#### Hybridization Buffer:

*Final composition is 10% (vol./vol.) formamide in Hybridization Buffer*

**Hybridization Buffer should be mixed fresh for each experiment:**

Due to viscosity of the solution, we recommend accounting for a 10% final volume excess in order to have enough Hybridization Buffer for all of your samples.

**For a final volume of 1 mL, mix:**

- 900 µL Stellaris RNA FISH Hybridization Buffer (Biosearch Technologies Cat# SMF-HB1-10)
- 100 µL Deionized Formamide

**NOTE:** Do not freeze Hybridization Buffer.

**WARNING!** Formamide is a teratogen that is easily absorbed through the skin and should be used in a chemical fume hood.

**WARNING!** Be sure to let the formamide warm to room temperature before opening the bottle.

#### Wash Buffer A (10 mL):

*Final composition is 10% (vol./vol.) formamide in 1X Wash Buffer A*

**Mix and dilute Wash Buffer A fresh for each experiment:**

**For a final volume of 10 mL, mix:**

- 2 mL Stellaris RNA FISH Wash Buffer A (Biosearch Technologies Cat# SMF-WA1-60)
- Add 7 mL Nuclease-free water
- Add 1 mL Deionized Formamide
- Mix well by vortexing gently.

#### Wash Buffer B:

**Add Nuclease-free water to Wash Buffer B bottle upon first use.**

Add 88 mL of Nuclease-free water to bottle (Biosearch Technologies Cat# SMF-WB1-20) before use. Mix thoroughly.

### **Nuclear Stain for use after hybridization:**

4',6-diamidino-2-phenylindole (DAPI) prepared in Wash Buffer A (see above) at 5 ng/mL. This solution is to be used in Step J below.

### **Mounting media:**

Vectashield Mounting Medium from Vector Laboratories (#H-1000).

**NOTE:** For best results, samples mounted with Vectashield Mounting Medium should be imaged the same day.

## **Protocol for *D. Melanogaster* Wing Imaginal Discs**

**NOTE:** This protocol is adapted from the Raj lab protocol<sup>3</sup> and has not been tested at Biosearch Technologies. Currently, we can only offer limited support for the use of Stellaris probes on *D. melanogaster*.

### **Fixation of *D. Melanogaster* Wing Imaginal Discs using chambered coverglass**

- a) Submerge 3rd instar larvae in 1 mL of 1X PBS and dissect to release wing imaginal discs.
- b) Place discs at the bottom of a chambered coverglass. The discs should readily stick.
- c) Gently aspirate 1X PBS, and then add 1mL of fixation buffer.
- d) Incubate at room temperature for 45 minutes.
- e) Aspirate fixation buffer, then wash twice with 1 mL of 1X PBS.
- f) To permeabilize, aspirate 1X PBS and add 1 mL of 70% ethanol and store overnight at +2 to +8 °C. Discs can be stored at +2 to +8 °C in 70% ethanol up to a week before hybridization.

### **Hybridization of *D. Melanogaster* Wing Imaginal Discs using chambered coverglass**

If frozen before using, warm the reconstituted probe solution to room temperature. Mix well by vortexing, then centrifuge briefly.

To prepare the Hybridization Buffer containing probe, add 1 µL of probe stock solution to 100 µL of Hybridization Buffer, and then vortex and centrifuge (enough for one coverglass). This creates a working probe solution of 125 nM. This solution will be used on step d.

- a) Aspirate the 70% ethanol off the chambered coverglass containing discs.
- b) Add 1 mL of Wash Buffer A (see recipe above), and incubate at room temperature for 2-5 minutes.
- c) Assemble humidified chamber: 150 mm tissue culture plate; place a single water-saturated paper towel alongside the inner chamber edge. This chamber will help prevent evaporation of the probe solution from under the coverglass.
- d) Aspirate Wash Buffer A, and dispense 100 µL of the Hybridization Buffer containing probe onto the discs within the chambered coverglass.
- e) Gently place a clean 18 x 18 mm #1 coverglass over the discs within the chambered coverglass.
- f) Place the chambered coverglass into the humidified chamber. Cover the humidified chamber with the lid and seal with Parafilm®.
- g) Incubate in the dark at 37 °C for at least 4 hours. (Incubation can be continued up to 16 hours).
- h) Add 1 mL of Wash Buffer A to the chambered coverglass. Use forceps to carefully remove the 18 x 18 mm coverglass so as not to disturb the discs underneath.
- i) Incubate in the dark at 37 °C for 30 minutes.
- j) Aspirate Wash Buffer A, and then add 1 mL of DAPI nuclear stain (Wash Buffer A consisting of 5 ng/mL DAPI) to counterstain the nuclei.
- k) Incubate in the dark at 37 °C for 30 minutes.
- l) Aspirate the DAPI staining buffer, and then add 1 mL of Wash Buffer B. Incubate at room temperature for 2-5 minutes.
- m) Add a small drop (approximately 15-30 µL) of Vectashield Mounting Medium onto the discs within the chambered coverglass. Place a clean 18 x 18 mm #1 coverglass over the discs to evenly spread the mounting medium.

### **Proceed to imaging.**

## **References**

1. Raj, A., van den Bogaard, P., Rifkin, S.A., van Oudenaarden, A., and Tyagi, S. Imaging individual mRNA molecules using multiple singly labeled probes. *Nat. Methods* 2008; 5, 877-9.
2. Femino, A.M., Fay, F.S., Fogarty, K., and Singer, R.H. Visualization of single RNA transcripts in situ. *Science* 1998; 280: 585-590.
3. Raj, A., Tyagi, S. Detection of individual endogenous RNA transcripts in situ using multiple singly labeled probes. *Methods in Enzymology*. 2010; 472: 365-86.

## Guidelines for Citing the Use of Stellaris RNA FISH Probes and Methods in Scientific Publications

Please acknowledge the use of Stellaris RNA FISH Probes and/or protocols in the experimental **Materials and Methods** or **Methods** section of your manuscript. Refer to the following examples as guidelines for proper citation of the Stellaris RNA FISH probe sets and/or protocols:

**Citing ShipReady and DesignReady Probe Sets:** "Stellaris® FISH Probes recognizing <catalogued gene set name> and labeled with <your dye of choice> (Catalog #, Biosearch Technologies, Inc., Petaluma, CA) were hybridized to <samples>, following the manufacturer's instructions available online at [www.biosearchtech.com/stellarisprotocols](http://www.biosearchtech.com/stellarisprotocols). Briefly, <describe any deviations from the published protocol or a short summary of what was actually performed>."

**Citing Custom Probe Sets Designed with the Stellaris® FISH Probe Designer:** "Custom Stellaris® FISH Probes were designed against <your RNA of interest (include NM# and nucleotides covered if relevant)> by utilizing the Stellaris® FISH Probe Designer (Biosearch Technologies, Inc., Petaluma, CA) available online at [www.biosearchtech.com/stellarisdesigner](http://www.biosearchtech.com/stellarisdesigner). The <samples> were hybridized with the <your RNA of interest> Stellaris FISH Probe set labeled with <your dye of choice> (Biosearch Technologies, Inc.), following the manufacturer's instructions available online at [www.biosearchtech.com/stellarisprotocols](http://www.biosearchtech.com/stellarisprotocols). Briefly, <describe any deviations from the published protocol or a short summary of what was actually performed>."

**Citing Custom Probe Sets Utilizing Previously Published Sequences:** "Custom Stellaris® FISH Probes recognizing <your RNA of interest (include NM# and nucleotides covered if relevant)> and labeled with <your dye of choice>, were purchased from Biosearch Technologies, Inc. (Petaluma, CA). Probe set sequences utilized in the experiments have been previously described <cite published manuscript>. The <samples> were hybridized with the <your RNA of interest> Stellaris® FISH Probe set, following the manufacturer's instructions available online at [www.biosearchtech.com/stellarisprotocols](http://www.biosearchtech.com/stellarisprotocols). Briefly, <describe any deviations from the published protocol or a short summary of what was actually performed>."

**Citing 3' Amine Oligos in Plates Used for Stellaris RNA FISH Designed with the Stellaris® FISH Probe Designer:** "Custom 3' amine oligos in plates were designed against <your RNA of interest (include NM# and nucleotides covered if relevant)> by utilizing the Stellaris® FISH Probe Designer (Biosearch Technologies, Inc., Petaluma, CA) available online at [www.biosearchtech.com/stellarisdesigner](http://www.biosearchtech.com/stellarisdesigner). Probes were labeled with <your dye of choice> using <insert your labeling protocol or citation of previously published labeling protocol>. The <samples> were hybridized with the <your RNA of interest> oligonucleotides (Biosearch Technologies, Inc.), following the manufacturer's instructions available online at [www.biosearchtech.com/stellarisprotocols](http://www.biosearchtech.com/stellarisprotocols). Briefly, <describe any deviations from the published protocol or a short summary of what was actually performed>."

**Citing 3' Amine Oligos in Plates Used for Stellaris RNA FISH Using Previously Published Sequences:** "Custom 3' amine oligos in plates recognizing <your RNA of interest (include NM# and nucleotides covered if relevant)> were purchased from Biosearch Technologies, Inc. (Petaluma, CA). Probe set sequences utilized in the experiments have been previously described <cite published manuscript>. Probes were labeled with <your dye of choice> using <insert your labeling protocol or citation of previously published labeling protocol>. The <samples> were hybridized with the <your RNA of interest> oligonucleotides (Biosearch Technologies, Inc.), following the manufacturer's instructions available online at [www.biosearchtech.com/stellarisprotocols](http://www.biosearchtech.com/stellarisprotocols). Briefly, <describe any deviations from the published protocol or a short summary of what was actually performed>."

## Technical Support

If you require additional information or technical assistance please feel free to e-mail our Technical Support Group at: [techsupport@biosearchtech.com](mailto:techsupport@biosearchtech.com). Our knowledgeable staff is also available for telephone consultation from 8:00 AM to 5:00 PM, Monday through Friday, Pacific Time. Please contact us at:

1.800.GENOME.1 (436.6631) US & Canada only

+1.415.883.8400 telephone

+1.415.883.8488 fax

2199 S. McDowell Blvd, Petaluma, CA 94954-6904 USA | [info@biosearchtech.com](mailto:info@biosearchtech.com)

©2015 Biosearch Technologies, Inc. All rights reserved.

This product is sold under license from Rutgers, The State University of New Jersey, and may be used under its patent rights for Research Use Only. Information on licensing programs for the commercial and/or diagnostic use of these products is available by contacting [licensing@biosearchtech.com](mailto:licensing@biosearchtech.com).

Biosearch, the Biosearch Technologies logo, and Stellaris are registered trademark of Biosearch Technologies, Inc., Products and technologies appearing in this brochure may have trademark or patent restrictions associated with them. Please see [www.biosearchtech.com/legal](http://www.biosearchtech.com/legal) for full legal disclosure.

UI-207269 Rev. 1.0  
Eff. Date: 13 Apr 2015