
Development of Ru(bpy)₃ CPG for Solid-phase Synthesis of 3'-labeled Oligonucleotides

Mary K. Johansson, Matthew H. Lyttle, Ron M. Cook



Objectives

The goal of this study is to develop a CPG solid support for DNA synthesis that is labeled with a ruthenium tris 2,2'-bipyridyl complex. We have prepared this CPG, made 3'-labeled oligonucleotides and performed some functional assays.

Introduction

Fluorescent dyes traditionally used as oligonucleotide labels (fluoresceins, rhodamines, cyanines) share some advantages(+) and disadvantages(-):

- + absorb light efficiently (ε > 50,000)
- + emit light efficiently ($\phi_{F} > 0.1$)
- +/- different degrees of photostability
- small Stokes shift (absorption and emission curves overlap)
- self-quenching
- must excite at long wavelengths to get near-IR emission


While Ru(bpy)₃ labels have the following characteristics:

- +/- moderate ε (ca. 15,000) and ϕ_F (ca. 0.1)
- + photostabl
- + large Stokes shift (excite at 460 nm and get emission at 610 nm)
- + not self-quenching
 - -> can use multiple labels
- + long fluorescence lifetimes (ca. 1 μs)
- + useful redox properties with applications including: electrochemistry, electrochemical detection, chemiluminescence

Properties of Ru Complexes


Structural

Ruthenium is a transition metal. Ru is below iron in the periodic table. Ru⁺² binds to six ligands to form an octahedral complex. These complexes are typically stable to organic synthesis conditions.




one enantiomer of Ru(bpy)₃

Common Bidentate Ligands:




Electronic

With organic dyes (FAM, Cy5, etc), the same excited state is involved in the absorption and emission of light:



The emission of light from ruthenium complexes is not simple fluorescence (i.e. more than singlet states are involved). When Ru(bpy)₃⁺² absorbs light at 450 nm, **Metal to Ligand Charge Transfer (MLCT)** occurs:

 $Ru^{+2}(bpy)_3 \longrightarrow Ru^{+3}(bpy)_2(bpy^{-1})$




Synthesis of CPG with a Ru Complex

Background

There have been several reports of labeling oligos with transition metals, most often as a 5'-modification. Grinstaff and coworkers at Duke University have prepared a Ru(bpy)₃ phosphoramidite. To our best knowledge, a CPG has not previously been prepared.

Preparation and Evaluation of 3'-Ru(bpy)₃ labeled Oligos

- The Ru complexes are stable to standard DNA synthesis conditions, including ammonia deprotection.
- * Obtained high yields of labeled oligos
- * Preparative reversed-phase (RP) HPLC purification gave >98% pure samples (according to both analytical anion-exchange and RP HPLC).
- * Oligos 3'-labeled with Ru complexes are photostable. They are also stable and luminesce at high temperatures.



A reversed-phase HPLC chromatogram of crude 3'Ru(bpy)₂(mcbpy)-β-actin (a 26-mer). The Ru(bpy)₃ CPG underwent standard DNA synthesis conditions and deprotection of 2 hours at 60 °C in concentrated ammonia.

A₂₆₀
0.7

B-actin-5'FAM
with 3T complement

0.6

0.7

B-actin-5'FAM
Temp (C)


Temp (C)


 $3'Ru(bpy)_2(mcbpy)-\beta$ -actin

with 3T complement

00

The melting temperatures for duplexes with the β-actin sequence (a 26-mer) are the same whether it is 3'-labeled with a ruthenium complex or 5'-labeled with FAM.





absorption and emission curves of the ruthenium complex are very well separated.

An endpoint assay of a dual-labeled probe, 5'-labeled with a Black Hole QuencherTM. The emission of the ruthenium complex depends on its distance from

the Black Hole

Quencher.

Conclusions

We have successfully prepared a Ru(bpy)₃-labeled CPG. This CPG has been used to make a variety of 3'-Ru(bpy)₃ labeled oligos in high yield. These labeled oligos are stable and have behaved well in complementation assays. It has also been shown, using a 5'quencher-3'reporter dual-labeled probe, that Ru(bpy)₃ efficiently transfers energy to a dark quencher, a Black Hole QuencherTM. The quenching of Ru(bpy)₃ by BHQ-2 probably occurs via the FRET mechanism.

References

Lakowicz, J. *Principles of Fluorescence*Spectroscopy; Plenum: New York, 1999;
Chapter 20.

Hu, X.; Smith, G. D.; Sykora, M.; Lee, S. J.; Grinstaff, M. W. *Inorg. Chem.* **2000**, 39, 2500, and reference contained therein.

For more information on BHQs or other Biosearch products, please visit booth 435 or our website www.biosearchtech.com.