

**REFERENCE DYES:**

**ROX dye (ROX):** provided with the Stratagene Brilliant Kit (standard).  
**Buffer A (BuA):** (ABI) this dye is in PCR buffer (standard).  
**3'ROX-Spacer 9 (RS9):** synthesized using ROX CPG & a PEG spacer amidite. The spacer is to increase solubility & stability.








**3'Cal Red-Spacer 9-Spacer 9 (CR18):** synthesized using Cal Red (Texas Red) CPG & two PEG spacer amidite couplings. The double spacer is to compare to the RS9 for increase in solubility & stability.  
**Blue 636 (B636):** from MegaBases. This dye is reported to emit from 630-640 nm.  
**5' 6FAM-T6-3' ROX (FT6R):** oligonucleotide synthesized to mimic BuA.

**Development and Evaluation of Passive Reference Candidate Compounds for Normalization of Signal in Real-Time PCR**


Cynthia Potter, Mary Katherine Johansson, Daren Dick, Ron Cook, Biosearch Technologies Inc., 81 Digital Dr., Novato CA, 94949-5750

**Purpose:** to test various reference dyes (RD) in real-time PCR.  
**Criterion:** >RD should tighten the threshold cycle (Ct) spread (threshold & baseline set by software).  
 >RD fluorescence intensity should remain constant in no template control (NTC), no amplification control (NAC) & all concentrations of template.  
 >Reporter (FAM) fluorescence intensity should remain constant in NTC & NAC

**Purpose for figures presented:**  
 > Multicomponent View (MV) is for viewing the fluorescence intensity of each dye over all PCR cycles  
 > Amplification plots are for calculation of Ct and standard deviation (SD) & viewing quality of amplification  
**Instrument platforms utilized:**  
 > ABI Prism 7700  
 > Corbett Research Rotor-Gene.  
**PCR conditions:**  
 400 nM ApoB primers, 200 nM ApoB FAM/BHQ-1 probe template: Human Genomic DNA (Clontech)  
**Initial run:** 100 ng, 10 ng, NTC  
 RD concentrations were adjusted so that their emission was less than the FAM emission.  
**Final run:** 100 ng, 10 ng, 1 ng, NTC



On the 7700, the RuC signals were constant, but the FAM signals increased steadily for the NAC & NTC samples. For amplification samples, the FAM signals increased steadily from the beginning of the PCR cycles and then increased with Ct, while the RuC signals remained constant. However, on the Rotor-Gene the RuC signals increased with Ct but remained quite constant for NAC and NTC samples. RuC1 and RuC2 performed similarly on the Rotor-Gene and RuC2 and RuC3 did not have such dramatic increases as RuC1, so RuC3 is expected to perform identically to RuC2 on the Rotor-Gene.



**Conclusions:**

The best reference dye was **ROX** for this study followed closely by **CR18**. This is based on Ct spread, signal reliability for many replicates of several concentrations, linearity of exponential portion of amplification curves.

**Buffer A** performed very well, but exhibited some ROX drop.

The **CR18** performed better than the **RS9**. This which could be attributed to the longer spacer on the **CR18**, which may prevent interaction with the PCR and improve stability. **RS9** performed very well--similar to **Buffer A**. A new compound consisting of **ROX** and two spacers may work as well as **CR18**.

**Blue 636**, with emission at ~645 nm, is still very interesting and warrants further study.

The difference in performance between the FRET oligos, **FT6R** and **Buffer A**, was unexpected and indicates that the dye in **Buffer A** is more stable for this particular assay. This merits further comparison, because this may be template or primer dependent.

The Ruthenium complexes, with their large Stokes' shift and long emission wavelength, show promise as reference dyes. However, more work is needed to overcome their increases in emission. Perhaps they are better suited as reporter dyes.

**Acknowledgements:**  
 We would like to thank Renee Horner for her assistance with the 7700 software. We would also like to thank Julie Vold for her assistance with running the 7700. For the use of the Rotor-Gene, we would like to thank Phenix Research Products, and especially Todd Deppe for his assistance with the software.