


Detection of Long Noncoding RNA in Cancer by Fluorescence *In Situ* Hybridization

Arturo V. Orjalo, Sally R. Coassin, Hans E. Johansson, and Jerry L. Ruth
81 Digital Drive, Novato, CA 94949-5728 USA

HYPOTHESIS

Long noncoding RNAs (lncRNAs) cannot be detected by immunohistochemistry, as no protein product is produced. Fluorescence *in situ* hybridization (FISH) using labeled DNA probes against chromosomal targets is also not useful, as transcription is not detected or quantified. The use of fluorescent synthetic oligonucleotide probes directed against the transcribed lncRNA should allow both the detection and quantification of the targets. To demonstrate this, we designed and tested fluorescent oligonucleotide probe sets against six lncRNAs (MALAT1, NEAT1, XIST, HOTAIR, H19, and PCA3) implicated in cancer.

METHODS



For each target, a mix of multiple 20mer oligonucleotides, each with a single Quasar® 570 fluorophore was designed* and synthesized. All probe sets contained ≥ 36 oligonucleotides each. Hybridizations were carried out on coverslips or slides for 4 to 16 hours at 37 °C in 100 μ L hybridization solution (10% dextran sulfate, 10% formamide in 2X SSC). Samples were then washed, nuclear stained with DAPI, and imaged.

* Stellaris® probe designer and protocols: www.biosearchtech.com/stellaris

RESULTS




Figure 1. NEAT1 and MALAT1 RNA FISH

(A) NEAT1 RNA distribution of its transcript variants in A549 cells. The short variant can be observed outside of nuclear paraspeckles, as seen in the long exposure.
(B) MALAT1 RNA distribution in nuclear speckles in A549 cells.




Figure 2. XIST RNA distribution in various human cell lines




Figure 3. HOTAIR and H19 RNA FISH

(A) HOTAIR RNA distribution in HeLa cells and HB2 cells.
(B) H19 RNA distribution and quantification in a single A549 cell.

REFERENCES

- Batish M, Raj A, Tyagi S. *Methods Mol Biol.* 2011;714:3-13.
- de Kok JB, Verhaegh GW, Roelofs RW, Hessels D, Kiemeny LA, Aalders TW, Swinkels DW, Schalken JA. *Cancer Res.* 2002 May 1;62(9):2695-8.
- Mao YS, Sunwoo H, Zhang B, Spector DL. *Nat Cell Biol.* 2011 Jan;13(1):95-101.
- Nakagawa S, Hirose T. *Cell Mol Life Sci.* 2012 Sep;69(18):3027-36.
- Nora EP, Lajoie BR, Schulz EG, Giorgetti L, Okamoto I, Servant N, Piolot T, van Berkum NL, Meisig J, Sedat J, Gribnau J, Barillot E, Blüthgen N, Dekker J, Heard E. *Nature.* 2012 Apr 11;485(7398):381-5.
- Rinn JL, Kertesz M, Wang JK, Squazzo SL, Xu X, Bruggmann SA, Goodnough LH, Helms JA, Farnham PJ, Segal E, Chang HY. *Cell.* 2007 Jun 29;129(7):1311-23.
- Wilusz JE, Freier SM, Spector DL. *Cell.* 2008 Nov 28;135(5):919-32.
- Yang F, Bi J, Xue X, Zheng L, Zhi K, Hua J, Fang G. *FEBS J.* 2012 Sep;279(17):3159-65.




Figure 4. Malat1 and Neat1 RNA FISH in FFPE mouse liver tissue
(A) Malat1 RNA probed in liver tissue knocked down for Malat1.
(B) Neat1 and beta-actin RNA distribution in liver tissue.




Figure 5. HOTAIR, H19, and PCA3 RNA FISH in human tissue
(A) HOTAIR and H19 RNA distribution in frozen ovarian tissue.
(B) PCA3 RNA distribution in FFPE prostate tissue (collected with appropriate permissions).

CONCLUSIONS

Long noncoding RNAs are easily detected by fluorescence *in situ* hybridization using a mixture of oligonucleotide probes, each labeled with a single fluorophore. Single molecule detection allows localization and quantification of the lncRNAs in question.

Stellaris® is a registered trademark of Biosearch Technologies, Inc. Products and technologies appearing in this poster may have trademark or patent restrictions associated with them. Please see www.biosearchtech.com/legal for full legal disclosure. This product is sold under license from the University of Medicine and Dentistry of New Jersey and may be used under its patent rights for Research Use Only.

Scan this QR Code to learn more about Stellaris RNA FISH Probes.

